

HEALTHCARE SHOULD FEEL HUMAN.

Holon Software Requirements

Introduction

The purpose of this document is to provide an overview of the software and minimum requirements to run on a workstation. Review the information below to ensure your organization's workstations meet the technical minimum requirements. If minimum requirements are met, please fill out the Practice Planner Form (link provided at the end of this document).

Holon Software

There are two pieces of software to support Holon Community: the Holon client and the Holon sensor.

Holon Client

The Holon client is a desktop application which enables end users to interact with features as scoped for the project.

The Holon client is to be installed on every user workstation that will require the functionality of the solution. **If admin rights are required for users to install apps on their computer, we will need IT involvement and participation during installation or any roll out planning.** For large organizations that support automated software deployment tools, Holon can provide a compatible version of the client installer. Instructions to install and log in will be provided during the project implementation.



Holon Sensor

The Holon sensor runs behind the scenes to detect the necessary information to allow the Holon client to work. Holon sensors enable Holon to glean user and patient context for use within the Holon client.

The Holon sensor installation will vary depending on EMR specifics. Completion of the Practice Planner Form is critical in allowing the implementation team to select the appropriate sensor.

Machine Requirements

Windows OS

The Holon client may be installed on virtually any recent Windows computer. The following information outlines recommended machine specifications:

- Windows 10 or higher (64 bit)
- Processor: 2 GHz or greater, a minimum CPU of i5 quad core or equivalent
- RAM: 8 GB min installed
- Disk space: 550MB

The Holon client may also be installed on thin client servers if needed.

Mac OS

Currently, there is no Mac support.

Screen Resolution

Screen resolution settings for workstations where the Holon client will be installed must have a minimum display resolution of 720p and 100% scaling.

Practice Planner Form

Link: <https://www.holonsolutions.com/practice-planner-form/>

Password: HolonPlanner2025

User Account Template:

<https://www.holonsolutions.com/wp-content/uploads/2025/02/UploadUserAccountsTemplate.xlsx>

